




VP39 Antibody

Product Code	CSB-PA890502XA01FRX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9T1Q9
Immunogen	Recombinant Acyrthosiphon pisum secondary endosymbiont phage 1 (Bacteriophage APSE-1) VP39 protein
Raised In	Rabbit
Species Reactivity	Acyrthosiphon pisum secondary endosymbiont phage 1 (Bacteriophage APSE-1)
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein A/G
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Acyrthosiphon pisum secondary endosymbiont phage 1 (Bacteriophage APSE-1)
Target Names	39

Image



Western Blot

Positive WB detected in Recombinant protein
 All lanes: VP39 antibody at 1:2000
 Secondary
 Goat polyclonal to rabbit IgG at 1/50000 dilution
 Predicted band size: 35.5 kDa
 Observed band size: 36 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.