



virD2 Antibody

Product Code	CSB-PA361912XA01RKV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06668
Immunogen	Recombinant Rhizobium radiobacter (Agrobacterium tumefaciens) (Agrobacterium radiobacter) virD2 protein
Raised In	Rabbit
Species Reactivity	Rhizobium radiobacter (Agrobacterium tumefaciens) (Agrobacterium radiobacter)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Rhizobium radiobacter (Agrobacterium tumefaciens) (Agrobacterium radiobacter)
Lead Time	Made-to-order (14-16 weeks)
Target Names	virD2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.