



old Antibody

Product Code	CSB-PA320238XA01FSA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13520
Immunogen	Recombinant Escherichia phage P2 (Bacteriophage P2) old protein
Raised In	Rabbit
Species Reactivity	Escherichia phage P2 (Bacteriophage P2)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Escherichia phage P2 (Bacteriophage P2)
Lead Time	Made-to-order (14-16 weeks)
Target Names	old
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.