



Intermediate capsid protein VP6 Antibody

Product Code	CSB-PA361082XA01FXW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03531
Immunogen	Recombinant Rotavirus A (strain RVA/SA11-Both/G3P5B[2]) (RV-A) (Simian Agent 11 (strain Both)) Intermediate capsid protein VP6 protein
Raised In	Rabbit
Species Reactivity	Rotavirus A (strain RVA/SA11-Both/G3P5B[2]) (RV-A) (Simian Agent 11 (strain Both))
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Rotavirus A (strain RVA/SA11-Both/G3P5B[2]) (RV-A) (Simian Agent 11 (strain Both))
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.