



xynA Antibody

Product Code	CSB-PA334790XA01CAL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P40944
Immunogen	Recombinant Caldicellulosiruptor sp. (strain Rt8B.4) xynA protein
Raised In	Rabbit
Species Reactivity	Caldicellulosiruptor sp. (strain Rt8B.4)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Caldicellulosiruptor sp. (strain Rt8B.4)
Lead Time	Made-to-order (14-16 weeks)
Target Names	xynA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.