



xInB Antibody

Product Code	CSB-PA345063XA01DOZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P55330
Immunogen	Recombinant Aspergillus niger xInB protein
Raised In	Rabbit
Species Reactivity	Aspergillus niger
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Aspergillus niger
Lead Time	Made-to-order (14-16 weeks)
Target Names	xInB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.