



y00F Antibody

Product Code	CSB-PA334559XA01EDZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P39416
Immunogen	Recombinant Enterobacteria phage T4 (Bacteriophage T4) y00F protein
Raised In	Rabbit
Species Reactivity	Enterobacteria phage T4 (Bacteriophage T4)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Enterobacteria phage T4 (Bacteriophage T4)
Lead Time	Made-to-order (14-16 weeks)
Target Names	y00F
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.