



frd.2 Antibody

Product Code	CSB-PA751875XA01FXX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q76VH1
Immunogen	Recombinant Enterobacteria phage LZ1 (Bacteriophage LZ1) frd.2 protein
Raised In	Rabbit
Species Reactivity	Enterobacteria phage LZ1 (Bacteriophage LZ1)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Enterobacteria phage LZ1 (Bacteriophage LZ1)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	frd.2