



XXXIII Antibody

Product Code	CSB-PA661000XA01EDO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q3T4N8
Immunogen	Recombinant Enterobacteria phage PRD1 (Bacteriophage PRD1) XXXIII protein
Raised In	Rabbit
Species Reactivity	Enterobacteria phage PRD1 (Bacteriophage PRD1)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Enterobacteria phage PRD1 (Bacteriophage PRD1)
Lead Time	Made-to-order (14-16 weeks)
Target Names	XXXIII
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.