



XYL1 Antibody

Product Code	CSB-PA845184XA01CZG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8X195
Immunogen	Recombinant Candida boidinii (Yeast) XYL1 protein
Raised In	Rabbit
Species Reactivity	Candida boidinii (Yeast)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Candida boidinii (Yeast)
Lead Time	Made-to-order (14-16 weeks)
Target Names	XYL1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.