



xynA Antibody

Product Code	CSB-PA334985XA01FUI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P23556
Immunogen	Recombinant Caldicellulosiruptor saccharolyticus (Caldocellum saccharolyticum) xynA protein
Raised In	Rabbit
Species Reactivity	Caldicellulosiruptor saccharolyticus (Caldocellum saccharolyticum)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Caldicellulosiruptor saccharolyticus (Caldocellum saccharolyticum)
Lead Time	Made-to-order (14-16 weeks)
Target Names	xynA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.