



frd.3 Antibody

Product Code	CSB-PA755258XA01FXG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q76VG7
Immunogen	Recombinant Enterobacteria phage LZ5 (Bacteriophage LZ5) frd.3 protein
Raised In	Rabbit
Species Reactivity	Enterobacteria phage LZ5 (Bacteriophage LZ5)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Enterobacteria phage LZ5 (Bacteriophage LZ5)
Lead Time	Made-to-order (14-16 weeks)
Target Names	frd.3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.