



Uncharacterized 5.8 kDa Antibody

Product Code	CSB-PA657801XA01FWP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q37928
Immunogen	Recombinant Escherichia phage HK022 (Bacteriophage HK022) Uncharacterized 5.8 kDa protein
Raised In	Rabbit
Species Reactivity	Escherichia phage HK022 (Bacteriophage HK022)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Escherichia phage HK022 (Bacteriophage HK022)
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.