



ZP2 Antibody

Product Code	CSB-PA027120XA01MOY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O77726
Immunogen	Recombinant Macaca radiata (Bonnet macaque) ZP2 protein
Raised In	Rabbit
Species Reactivity	Macaca radiata (Bonnet macaque)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Macaca radiata (Bonnet macaque)
Lead Time	Made-to-order (14-16 weeks)
Target Names	ZP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.