



Uncharacterized 5.5 kDa protein in replication origin region Antibody

Product Code	CSB-PA321256XA01ENL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14505
Immunogen	Recombinant Escherichia coli Uncharacterized 5.5 kDa protein in replication origin region protein
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Escherichia coli
Lead Time	Made-to-order (14-16 weeks)