



Uncharacterized protein in scsB 5'region Antibody

Product Code	CSB-PA188148XA01TNR
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25124
Immunogen	Recombinant Thermus thermophilus Uncharacterized protein in scsB 5'region protein
Raised In	Rabbit
Species Reactivity	Thermus thermophilus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Thermus thermophilus
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.