



zraS Antibody

Product Code	CSB-PA860569XA01KBE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9APE0
Immunogen	Recombinant Klebsiella oxytoca zraS protein
Raised In	Rabbit
Species Reactivity	Klebsiella oxytoca
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Klebsiella oxytoca
Lead Time	Made-to-order (14-16 weeks)
Target Names	zraS
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.