



## yhfG Antibody

|                     |                                                                                       |
|---------------------|---------------------------------------------------------------------------------------|
| Product Code        | CSB-PA334318XA01SXB                                                                   |
| Storage             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                         |
| Uniprot No.         | P37771                                                                                |
| Immunogen           | Recombinant Salmonella typhimurium (strain LT2 / SGSC1412 / ATCC 700720) yhfG protein |
| Raised In           | Rabbit                                                                                |
| Species Reactivity  | Salmonella typhimurium (strain LT2 / SGSC1412 / ATCC 700720)                          |
| Tested Applications | ELISA, WB (ensure identification of antigen)                                          |
| Form                | Liquid                                                                                |
| Conjugate           | Non-conjugated                                                                        |
| Storage Buffer      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4      |
| Purification Method | Antigen Affinity Purified                                                             |
| Isotype             | IgG                                                                                   |
| Clonality           | Polyclonal                                                                            |
| Product Type        | Polyclonal Antibody                                                                   |
| Immunogen Species   | Salmonella typhimurium (strain LT2 / SGSC1412 / ATCC 700720)                          |
| Lead Time           | Made-to-order (14-16 weeks)                                                           |
| Gene Names          | yhfG                                                                                  |