



Uncharacterized 65 kDa protein in hyaluronidase region Antibody

Product Code	CSB-PA324458XA01SMS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15317
Immunogen	Recombinant Streptococcus pyogenes phage H4489A Uncharacterized 65 kDa protein in hyaluronidase region protein
Raised In	Rabbit
Species Reactivity	Streptococcus pyogenes phage H4489A
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Streptococcus pyogenes phage H4489A
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.