



Uncharacterized 3.7 kDa protein in rusA 5'region Antibody

Product Code	CSB-PA655011XA01ECB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q37869
Immunogen	Recombinant Enterobacteria phage 82 (Bacteriophage 82) Uncharacterized 3.7 kDa protein in rusA 5'region protein
Raised In	Rabbit
Species Reactivity	Enterobacteria phage 82 (Bacteriophage 82)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Enterobacteria phage 82 (Bacteriophage 82)
Lead Time	Made-to-order (14-16 weeks)