



Uncharacterized metal-dependent hydrolase in metS 3'region Antibody

Product Code	CSB-PA560311XA01GFM
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P23918
Immunogen	Recombinant Geobacillus stearothermophilus (Bacillus stearothermophilus) Uncharacterized metal-dependent hydrolase in metS 3'region protein
Raised In	Rabbit
Species Reactivity	Geobacillus stearothermophilus (Bacillus stearothermophilus)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Geobacillus stearothermophilus (Bacillus stearothermophilus)
Lead Time	Made-to-order (14-16 weeks)