



AC4 Antibody

Product Code	CSB-PA315128XA01AAE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0C6F9
Immunogen	Recombinant Abutilon mosaic virus (isolate West India) (AbMV) AC4 protein
Raised In	Rabbit
Species Reactivity	Abutilon mosaic virus (isolate West India) (AbMV)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Abutilon mosaic virus (isolate West India) (AbMV)
Lead Time	Made-to-order (14-16 weeks)
Target Names	AC4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.