



FATB1 Antibody

Product Code	CSB-PA886743XA01GHB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9SQI3
Immunogen	Recombinant Gossypium hirsutum (Upland cotton) (Gossypium mexicanum) FATB1 protein
Raised In	Rabbit
Species Reactivity	Gossypium hirsutum (Upland cotton) (Gossypium mexicanum)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Gossypium hirsutum (Upland cotton) (Gossypium mexicanum)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	FATB1