



HB2 Antibody

Product Code	CSB-PA842417XA01GHB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q93Y92
Immunogen	Recombinant Gossypium hirsutum (Upland cotton) (Gossypium mexicanum) HB2 protein
Raised In	Rabbit
Species Reactivity	Gossypium hirsutum (Upland cotton) (Gossypium mexicanum)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Gossypium hirsutum (Upland cotton) (Gossypium mexicanum)
Lead Time	Made-to-order (14-16 weeks)
Target Names	HB2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.