



Probable phospholipid hydroperoxide glutathione peroxidase Antibody

Product Code	CSB-PA526179XA01GHB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O49069
Immunogen	Recombinant Gossypium hirsutum (Upland cotton) (Gossypium mexicanum) Probable phospholipid hydroperoxide glutathione peroxidase protein
Raised In	Rabbit
Species Reactivity	Gossypium hirsutum (Upland cotton) (Gossypium mexicanum)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Gossypium hirsutum (Upland cotton) (Gossypium mexicanum)
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.