



matK Antibody

Product Code	CSB-PA822073XA01FX Y
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q95EF3
Immunogen	Recombinant Kokia drynarioides (Hawaiian tree cotton) (Gossypium drynarioides) matK protein
Raised In	Rabbit
Species Reactivity	Kokia drynarioides (Hawaiian tree cotton) (Gossypium drynarioides)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Kokia drynarioides (Hawaiian tree cotton) (Gossypium drynarioides)
Lead Time	Made-to-order (14-16 weeks)
Target Names	matK
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.