



BC1 Antibody

Product Code	CSB-PA325182XA01AAE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21946
Immunogen	Recombinant Abutilon mosaic virus (isolate West India) (AbMV) BC1 protein
Raised In	Rabbit
Species Reactivity	Abutilon mosaic virus (isolate West India) (AbMV)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Abutilon mosaic virus (isolate West India) (AbMV)
Lead Time	Made-to-order (14-16 weeks)
Target Names	BC1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.