




Non-specific lipid-transfer Antibody

Product Code	CSB-PA670383XA01GHB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q43129
Immunogen	Recombinant Gossypium hirsutum (Upland cotton) (Gossypium mexicanum) Non-specific lipid-transfer protein (27-120aa)
Raised In	Rabbit
Species Reactivity	Gossypium hirsutum (Upland cotton) (Gossypium mexicanum)
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Gossypium hirsutum (Upland cotton) (Gossypium mexicanum)
Target Names	Non-specific lipid-transfer Antibody

Image



Western Blot

Positive WB detected in: recombinant protein
 All lanes: Non-specific lipid-transfer Antibody at 1:1000
 Secondary
 Goat polyclonal to rabbit IgG at 1/50000 dilution
 Predicted band size: 24.6 kDa
 Observed band size: 28 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.