



psbJ Antibody

Product Code	CSB-PA371608XA01GEP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A0ZZ49
Immunogen	Recombinant Gossypium barbadense (Sea-island cotton) (Egyptian cotton) psbJ protein
Raised In	Rabbit
Species Reactivity	Gossypium barbadense (Sea-island cotton) (Egyptian cotton)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Gossypium barbadense (Sea-island cotton) (Egyptian cotton)
Lead Time	Made-to-order (14-16 weeks)
Target Names	psbJ
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.