



ORF2a Antibody

Product Code	CSB-PA308580XA01TFN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P89678
Immunogen	Recombinant Tobacco streak virus (strain WC) (TSV) ORF2a protein
Raised In	Rabbit
Species Reactivity	Tobacco streak virus (strain WC) (TSV)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Tobacco streak virus (strain WC) (TSV)
Lead Time	Made-to-order (14-16 weeks)
Target Names	ORF2a
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.