



TIC214 Antibody

Product Code	CSB-PA645692XA01GHB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q2L951
Immunogen	Recombinant Gossypium hirsutum (Upland cotton) (Gossypium mexicanum) TIC214 protein
Raised In	Rabbit
Species Reactivity	Gossypium hirsutum (Upland cotton) (Gossypium mexicanum)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Gossypium hirsutum (Upland cotton) (Gossypium mexicanum)
Lead Time	Made-to-order (14-16 weeks)
Target Names	TIC214
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.