



Arginine biosynthesis bifunctional protein ArgJ, chloroplastic Antibody

Product Code	CSB-PA499543XA01ZAX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	C0PF72
Immunogen	Recombinant Zea mays (Maize) Arginine biosynthesis bifunctional protein ArgJ, chloroplastic protein
Raised In	Rabbit
Species Reactivity	Zea mays (Maize)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Zea mays (Maize)
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.