




Endochitinase B Antibody

Product Code	CSB-PA333407XA01ZAX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29023
Immunogen	Recombinant Zea mays (Maize) CHIB protein
Raised In	Rabbit
Species Reactivity	Zea mays (Maize)
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein A/G
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Zea mays (Maize)

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: Endochitinase B antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 33.8 kDa

Observed band size: 33 kDa