



GAPC2 Antibody

Product Code	CSB-PA276518XA01ZAX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q09054
Immunogen	Recombinant Zea mays (Maize) GAPC2 protein
Raised In	Rabbit
Species Reactivity	Zea mays (Maize)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Zea mays (Maize)
Lead Time	Made-to-order (14-16 weeks)
Target Names	GAPC2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.