



GRMZM2G118515 Antibody

Product Code	CSB-PA564440XA01ZAX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B4FIJ0
Immunogen	Recombinant Zea mays (Maize) GRMZM2G118515 protein
Raised In	Rabbit
Species Reactivity	Zea mays (Maize)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Zea mays (Maize)
Lead Time	Made-to-order (14-16 weeks)
Target Names	GRMZM2G118515
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.