



ORF3 Antibody

Product Code	CSB-PA362708XA01BER
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09510
Immunogen	Recombinant Barley yellow dwarf virus (isolate PAV) (BYDV) ORF3 protein
Raised In	Rabbit
Species Reactivity	Barley yellow dwarf virus (isolate PAV) (BYDV)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Barley yellow dwarf virus (isolate PAV) (BYDV)
Lead Time	Made-to-order (14-16 weeks)
Target Names	ORF3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.