



AR1 Antibody

Product Code	CSB-PA896663XA01FYE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9YPS5
Immunogen	Recombinant Mungbean yellow mosaic virus (strain Vigna) (MYMV) AR1 protein
Raised In	Rabbit
Species Reactivity	Mungbean yellow mosaic virus (strain Vigna) (MYMV)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Mungbean yellow mosaic virus (strain Vigna) (MYMV)
Lead Time	Made-to-order (14-16 weeks)
Target Names	AR1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.