



Os06g0520600 Antibody

Product Code	CSB-PA726464XA01OFG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5Z5Q3
Immunogen	Recombinant Oryza sativa subsp. japonica (Rice) Os06g0520600 protein
Raised In	Rabbit
Species Reactivity	Oryza sativa subsp. japonica (Rice)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Oryza sativa subsp. japonica (Rice)
Lead Time	Made-to-order (14-16 weeks)
Target Names	Os06g0520600
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.