



# EXPA22 Antibody

<b>Product Code</b>	CSB-PA686897XA01OFG
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q4PR44
<b>Immunogen</b>	Recombinant Oryza sativa subsp. japonica (Rice) EXPA22 protein
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Oryza sativa subsp. japonica (Rice)
<b>Tested Applications</b>	ELISA, WB (ensure identification of antigen)
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Oryza sativa subsp. japonica (Rice)
<b>Lead Time</b>	Made-to-order (14-16 weeks)
<b>Gene Names</b>	EXPA22