



Os02g0491800 Antibody

Product Code	CSB-PA753745XA01OFG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6K5P8
Immunogen	Recombinant Oryza sativa subsp. japonica (Rice) Os02g0491800 protein
Raised In	Rabbit
Species Reactivity	Oryza sativa subsp. japonica (Rice)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Oryza sativa subsp. japonica (Rice)
Lead Time	Made-to-order (14-16 weeks)
Target Names	Os02g0491800
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.