



Os05g0361200 Antibody

Product Code	CSB-PA607733XA01OFG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q0DIV0
Immunogen	Recombinant Oryza sativa subsp. japonica (Rice) Os05g0361200 protein
Raised In	Rabbit
Species Reactivity	Oryza sativa subsp. japonica (Rice)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Oryza sativa subsp. japonica (Rice)
Lead Time	Made-to-order (14-16 weeks)
Target Names	Os05g0361200
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.