



Outer capsid protein P8 Antibody

Product Code	CSB-PA335434XA01RMJ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29077
Immunogen	Recombinant Rice gall dwarf virus (RGDV) Outer capsid protein P8 protein
Raised In	Rabbit
Species Reactivity	Rice gall dwarf virus (RGDV)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Rice gall dwarf virus (RGDV)
Lead Time	Made-to-order (14-16 weeks)