



Outer capsid protein P3 Antibody

Product Code	CSB-PA847457XA01FWE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q98630
Immunogen	Recombinant Rice dwarf virus (isolate Fujian) (RDV) Outer capsid protein P3 protein
Raised In	Rabbit
Species Reactivity	Rice dwarf virus (isolate Fujian) (RDV)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Rice dwarf virus (isolate Fujian) (RDV)
Lead Time	Made-to-order (14-16 weeks)