



Minor outer capsid protein P2 Antibody

Product Code	CSB-PA529995XA01RMI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O55519
Immunogen	Recombinant Rice dwarf virus (isolate O) (RDV) Minor outer capsid protein P2 protein
Raised In	Rabbit
Species Reactivity	Rice dwarf virus (isolate O) (RDV)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Rice dwarf virus (isolate O) (RDV)
Lead Time	Made-to-order (14-16 weeks)