



Protein P1 Antibody

Product Code	CSB-PA329730XA01RER
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P27500
Immunogen	Recombinant Rice tungro bacilliform virus (isolate Philippines) (RTBV) Protein P1 protein
Raised In	Rabbit
Species Reactivity	Rice tungro bacilliform virus (isolate Philippines) (RTBV)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Rice tungro bacilliform virus (isolate Philippines) (RTBV)
Lead Time	Made-to-order (14-16 weeks)