



Minor outer capsid protein P9 Antibody

Product Code	CSB-PA801061XA01FWE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q85440
Immunogen	Recombinant Rice dwarf virus (isolate Fujian) (RDV) Minor outer capsid protein P9 protein
Raised In	Rabbit
Species Reactivity	Rice dwarf virus (isolate Fujian) (RDV)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Rice dwarf virus (isolate Fujian) (RDV)
Lead Time	Made-to-order (14-16 weeks)