



KIN5A Antibody

Product Code	CSB-PA906383XA01OFG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5W7C6
Immunogen	Recombinant Oryza sativa subsp. japonica (Rice) KIN5A protein
Raised In	Rabbit
Species Reactivity	Oryza sativa subsp. japonica (Rice)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Oryza sativa subsp. japonica (Rice)
Lead Time	Made-to-order (14-16 weeks)
Target Names	KIN5A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.