



Os10g0478200 Antibody

Product Code	CSB-PA768678XA01OFG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7XDC8
Immunogen	Recombinant Oryza sativa subsp. japonica (Rice) Os10g0478200 protein
Raised In	Rabbit
Species Reactivity	Oryza sativa subsp. japonica (Rice)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Oryza sativa subsp. japonica (Rice)
Lead Time	Made-to-order (14-16 weeks)
Target Names	Os10g0478200
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.