



ORF4 Antibody

Product Code	CSB-PA357694XA01BER
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09513
Immunogen	Recombinant Barley yellow dwarf virus (isolate PAV) (BYDV) ORF4 protein
Raised In	Rabbit
Species Reactivity	Barley yellow dwarf virus (isolate PAV) (BYDV)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Barley yellow dwarf virus (isolate PAV) (BYDV)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	ORF4