



Protein P7 Antibody

Product Code	CSB-PA322167XA01RMI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P22473
Immunogen	Recombinant Rice dwarf virus (isolate O) (RDV) Protein P7 protein
Raised In	Rabbit
Species Reactivity	Rice dwarf virus (isolate O) (RDV)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Rice dwarf virus (isolate O) (RDV)
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.